

Datasheet: 0100-0035

Description:	NATIVE LIPOTEICHOIC ACID
Name:	LIPOTEICHOIC ACID
Other names:	LTA
Format:	Purified
Product Type:	Antigen
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			*

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Purified lipoteichoic acid - liquid
Preparation	Phenolic extraction
Preservative Stabilisers	None present.
Approx. Protein Concentrations	1.0 mg/ml

Product Information Native lipoteichoic acid (LTA) purified from *Staphylococcus aureus*.

LTA is a linear polymer of phosphodiester-linked glycerol phosphate covalently attached to lipid. The lipid binds hydrophobically to the cell membrane and the polyglycerol phosphate portion extends into the cell wall. The biological role of LTA is not fully understood but when expressed on the cell exterior it is thought to act as an adhesin.

CAS Number: 56411-57-5

ELISA	*When coating use 10 ng - 1µg/well.
Storage	Store at -20°C only. Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 2 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation available at:
Material Safety Datasheet Documentation #10078 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10078.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M313264:180319'

Printed on 20 Mar 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)