

Datasheet: 106005

Description:	GOAT ANTI MOUSE LAMBDA LIGHT CHAIN:HRP
Specificity:	LAMBDA LIGHT CHAIN
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/4000 - 1/8000
Western Blotting			▪	
Immunoblotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Ig Fraction conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation Antisera to mouse lambda light chains were raised by repeated immunisations of goats with highly purified antigen. Purified Ig was prepared from whole serum by affinity chromatography on pooled mouse Igs with lambda chains covalently linked to agarose.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	50% glycerol
---------------------------------	--------------

Immunogen	Mouse paraproteins with lambda 1 and lambda 2 light chains.
------------------	---

External Database Links

UniProt:
[P01843](#) [Related reagents](#)
[P01844](#) [Related reagents](#)

Entrez Gene:

[433053](#) LOC433053 [Related reagents](#)

[110786](#) Iglc2 [Related reagents](#)

Specificity **Goat anti Mouse Lambda Light Chain antibody** recognizes mouse lambda 1 and lambda 2 light chains (free and bound) as demonstrated by ELISA. To minimise cross reactivity this product has been cross absorbed against mouse paraproteins with kappa light chains.

References

1. Rops, A.L. *et al.* (2010) The tetraspanin CD37 protects against glomerular IgA deposition and renal pathology. [Am J Pathol. 176 \(5\): 2188-97.](#)
2. Lewis, A. *et al.* (2017) Recovery from desensitization of IgE-dependent responses in human lung mast cells. [Clin Exp Allergy. Feb 25. \[Epub ahead of print\]](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted.

Shelf Life Please see label for expiry date

Health And Safety Information Material Safety Datasheet documentation available at:
Material Safety Datasheet Documentation #10106 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10106.pdf>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M314765:180412'

Printed on 13 Apr 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)