

Datasheet: 12004157

Description:	GOAT ANTI-MOUSE IgG StarBright™ Blue 700
Specificity:	IgG
Format:	StarBright Blue 700
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	5 x 400µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/2,500 - 1/5,000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Mouse

Product Form

Purified IgG conjugated to StarBright Blue 700 - lyophilized

Reconstitution

Resuspend contents of the tube in the indicated volume of distilled or deionized water and leave on ice for at least 30 min prior to use. Brief centrifugation (pulse spin for 2-3 sec at max speed in a tabletop microcentrifuge) may be used to collect the contents at the bottom of the tube. Do not centrifuge more than 10 sec.

Antiserum Preparation Antiserum to Mouse IgG was raised by repeated immunisation of Goats with highly purified antigen.

Buffer Solution

Phosphate buffered saline

Preservative

>80% Sucrose

Stabilisers

<1% Polyethylene Glycol

<0.04% 2-Methyl-2H-Isothiazol-3-One

<10% Bovine Serum Albumin

Product Information

StarBright Blue 700 Secondary Antibodies are ideal for fluorescent western blotting — either for the detection of a single target protein or for multiplex detection of several proteins on one blot, without stripping and reprobing. The StarBright Fluorophore is composed of a condensed polymer made up of multiple light-absorbing and -emitting monomers which provides an exceptionally bright signal (high quantum yield) compared to most traditional fluorophores.

StarBright Blue 700 Fluorescent Secondary Antibodies can be used with traditional fluorophores

like the DyLight 800 Fluorophore for multiplexing. In addition, StarBright Antibodies can be used with Bio-Rad's stain-free technology and/or the hFAB™ Rhodamine Anti-Housekeeping Primary Antibodies for protein normalization. These antibodies are optimized for use with the ChemiDoc™ MP Imaging System, permitting detection of up to three proteins in a single blot.

Features and Benefits

Easy multiplexing — StarBright Secondary Antibodies can be used with IRDye 800 or other common fluorophores to detect multiple proteins on the same blot

Short exposure time — the exceptional brightness of the StarBright Fluorophore leads to shorter exposure times compared to traditional fluorescent dyes: seconds vs. minutes

High signal-to-noise ratio — StarBright Secondary Antibodies emit in the red/far red, where there is minimal background ($Ex_{max}/Em_{max} = 470\text{ nm}/700\text{ nm}$)

Low nonspecific binding — StarBright Fluorophore is conjugated to highly cross-adsorbed antibodies

Goat Anti-Mouse IgG StarBright Blue 700 has been cross-adsorbed against bovine, goat, human, rabbit and rat.

Western Blotting	StarBright Blue 700 Secondary Antibodies can be used for multiplex fluorescent western blotting with the right combination of fluorophores. These fluorescent secondary antibodies can also be used to detect post-translational modifications like phosphorylation of a target protein. Note: StarBright Antibodies are not suitable for stripping and reprobing.
Instructions For Use	Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/12004157.pdf
Storage	1 year at -20°C lyophilized; 6 months at +4°C after resuspension. This product is photosensitive and should be protected from light DO NOT FREEZE the solubilized material
Shelf Life	Please see label for expiry date
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet documentation #12004157 available at https://www.bio-rad-antibodies.com/uploads/MSDS/12004157.pdf
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[GOAT ANTI-RABBIT IgG StarBright™ Blue 700 \(12004160\)](#)
[GOAT ANTI-RABBIT IgG StarBright™ Blue 700 \(12004161\)](#)
[GOAT ANTI-RABBIT IgG StarBright™ Blue 700 \(12004162\)](#)
[ANTI-ACTIN hFAB™ RHODAMINE ANTIBODY \(12004163\)](#)
[ANTI-ACTIN hFAB™ RHODAMINE ANTIBODY \(12004164\)](#)
[ANTI-TUBULIN hFAB™ RHODAMINE ANTIBODY \(12004165\)](#)
[ANTI-TUBULIN hFAB™ RHODAMINE ANTIBODY \(12004166\)](#)
[ANTI-GAPDH hFAB™ RHODAMINE ANTIBODY \(12004167\)](#)
[ANTI-GAPDH hFAB™ RHODAMINE ANTIBODY \(12004168\)](#)
[SHEEP ANTI RABBIT IgG:DyLight®800 \(STAR36D800GA\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com
'M310399:170921'

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

Printed on 13 Jan 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)