

Datasheet: 205004

| | |
|----------------------|--------------------------------|
| Description: | GOAT ANTI HUMAN IgA:Alk. Phos. |
| Specificity: | IgA |
| Format: | Alk. Phos. |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | | | ▪ | |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | ▪ | | | 1/2000 - 1/4000 |
| Immunoprecipitation | | | ▪ | |
| Western Blotting | | | ▪ | |
| Immunoblotting | | | ▪ | |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

| | |
|------------------------------|--|
| Target Species | Human |
| Product Form | Purified IgG conjugated to Alkaline Phosphatase - liquid |
| Preparation | Purified IgG prepared by affinity chromatography |
| Antiserum Preparation | Antisera to human IgA were raised by repeated immunisations of goats with highly purified antigen. |
| Buffer Solution | TRIS buffered saline, 1mM MgCl ₂ |
| Preservative | 0.1% Sodium Azide |
| Stabilisers | 50% Glycerol |
| Immunogen | Human IgA. |

External Database Links

UniProt:

[P01876](#)

[Related reagents](#)

[P01877](#)

[Related reagents](#)

Entrez Gene:[3493](#) IGHA1 [Related reagents](#)[3494](#) IGHA2 [Related reagents](#)

Specificity

Goat anti Human IgA antibody recognizes the heavy chain of human IgA as demonstrated by ELISA and Flow cytometry. IgA is the major category of antibody present in mucosal secretions and acts as a first line of defense against inhaled or ingested pathogens at the mucosal interface. IgA is also present in serum where it acts to eliminate pathogens via interaction with FcαR, also known as CD89. ([Woof et al. 2004](#)).
Biotinylated Goat anti Human IgA antibody (**205008**) has been cross-adsorbed against human IgD, IgG and IgM to ensure class specificity.

References

1. Wright, A.K. *et al.* (2012) Human Nasal Challenge with *Streptococcus pneumoniae* Is Immunising in the Absence of Carriage. [PLoS Pathog. 8: e1002622](#).
2. Pennington, S.H. *et al.* (2016) Polysaccharide-specific Memory B-cells Predict Protection Against Experimental Human Pneumococcal Carriage. [Am J Respir Crit Care Med. Jul 12. \[Epub ahead of print\]](#)

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

Please see label for expiry date.

Health And Safety Information

Material Safety Datasheet documentation #10322 available at:
10322: <https://www.bio-rad-antibodies.com/uploads/MSDS/10322.pdf>

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M314072:180412'

Printed on 01 May 2018