

Datasheet: 4750-0254

Description:	SHEEP ANTI HUMAN GROWTH HORMONE
Specificity:	GROWTH HORMONE
Other names:	SOMATOTROPIN
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.25 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Radioimmunoassays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Serum - liquid
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Purified human growth hormone.
External Database Links	<p>UniProt: P01241 Related reagents</p> <p>Entrez Gene: 2688 GH1 Related reagents</p>

Specificity

Sheep anti Human growth hormone antibody recognizes human growth hormone. Growth hormone is a small peptide hormone released by the anterior pituitary under the control of the hypothalamus. It acts to control body growth and metabolic activity.

References

1. Rowe, M. *et al.* (1994) Upregulation of bcl-2 by the Epstein-Barr virus latent membrane protein LMP1: a B-cell-specific response that is delayed relative to NF-kappa B activation and to induction of cell surface markers. [J Virol. 68: 5602-12](#).

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 18 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation available at:
Material Safety Datasheet Documentation #10081 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10081.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

Donkey Anti Sheep IgG (STAR88...) [DyLight@488](#), [DyLight@549](#), [DyLight@649](#),
[FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M310138:170914'

Printed on 14 Sep 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)