

Datasheet: 4890-0004P

Description:	SHEEP ANTI HUMAN HAPTOGLOBIN:HRP
Specificity:	HAPTOGLOBIN
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1/100 - 1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species Human

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to Haptoglobin were raised by repeated immunisations of sheep with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.01% Thiomersal

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen Native, purified human haptoglobin from serum.

External Database Links

UniProt:

[P00738](#) [Related reagents](#)

Entrez Gene:

[3240](#) HP [Related reagents](#)

Specificity

Sheep anti Human Haptoglobin antibody detects haptoglobin, also known as Hp, an acute phase plasma glycoprotein synthesized by the liver in response to infection or inflammation.

Haptoglobin binds haemoglobin released by erythrocytes with high affinity, and functions to block the oxidative activity of haemoglobin and prevent renal damage. The haptoglobin-haemoglobin complex is rapidly removed from the circulation by monocytes and macrophages, via the receptor CD163.

The haptoglobin genotype has been demonstrated to be an independent risk factor for cardiovascular disease in multiple epidemiological studies.

References

1. Carter, K. & Worwood, M. (2007) Haptoglobin: a review of the major allele frequencies worldwide and their association with diseases. [Int J Lab Hematol. 29 \(2\): 92-110.](#)
2. Guetta, J. *et al.* (2007) Haptoglobin genotype modulates the balance of Th1/Th2 cytokines produced by macrophages exposed to free hemoglobin. [Atherosclerosis. 191: 48-53.](#)
3. Araki, H. *et al.* (2015) Haptoglobin promoter polymorphism rs5472 as a prognostic biomarker for peptide vaccine efficacy in castration-resistant prostate cancer patients. [Cancer Immunol Immunother. 64 \(12\): 1565-73.](#)
4. Soejima, M. *et al.* (2014) Genetic factors associated with serum haptoglobin level in a Japanese population. [Clin Chim Acta. 433: 54-7.](#)

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

Health And Safety Information

Material Safety Datasheet documentation #10094 available at:
10094: <https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf>

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M275824:151012'

Printed on 01 May 2018