

Datasheet: 7370-9902

Description:	NATIVE RABBIT PHOSPHORYLASE B
Name:	PHOSPHORYLASE B
Format:	Purified
Product Type:	Purified Protein
Quantity:	25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Purified enzyme from rabbit skeletal muscle - liquid
Preservative Stabilisers	None present
Approx. Protein Concentrations	5.0 mg/ml

Product Information **Native Rabbit Phosphorylase B** is a member of the protein kinase superfamily. Phosphorylases catalyse the breakdown of glycogen in liver and muscle. The phosphorylase of skeletal muscle occurs in two forms, the active form (phosphorylase-a) and a much less active form (phosphorylase b). Phosphorylase b has a molecular weight of approximately 185kDa.

Activity	55-60 U/mg protein One unit will liberate one μ M of glucose-1-phosphate from glycogen and orthophosphate per minute at 37°C and pH 6.8.
Purity	Protein: 95%
Further Reading	1. Fischer, E.H. and Krebs, E.G. (1958) The isolation and crystallization of rabbit skeletal muscle phosphorylase b. J Biol Chem. 231: 65-71.
Storage	Store at -20°C only.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 18 months from date of despatch.

Health And Safety Information Material Safety Datasheet Documentation #10078 available at: <https://www.bio-rad-antibodies.com/uploads/MSDS/10078.pdf>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[SHEEP ANTI RABBIT PHOSPHORYLASE B \(7370-9808\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M298763:161209'

Printed on 06 Jun 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)