

## Datasheet: 9220-1995

<b>Description:</b>	NATIVE PIG TRYPSIN
<b>Name:</b>	TRYPSIN
<b>Format:</b>	Semi Pure
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	500 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Pig
<b>Product Form</b>	Partially purified enzyme from porcine pancreas - liquid
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	Total protein concentration 100 mg/ml
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P00761</a> <a href="#">Related reagents</a>

**Product Information** 9220-1995 is partially purified porcine pancreatic trypsin, a serine protease that shows preferential cleavage at Arg-|-Xaa, and Lys-|-Xaa sites.

**Activity** 1500 U/mg solid

1 unit causes an increase in absorbance of 0.0003 per minute at 25°C and 253 nm, pH 7.6 on the substrate benzoyl-L arginine ester (BAEE).

**Storage** Store at -20°C only.  
 Storage in frost-free freezers is not recommended.  
 This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend

microcentrifugation before use.

---

**Shelf Life** 18 months from date of despatch.

---

**Health And Safety Information** Material Safety Datasheet Documentation #10078 available at:  
<https://www.bio-rad-antibodies.com/uploads/MSDS/10078.pdf>

---

**Regulatory** For research purposes only

---

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M298779:161209'

**Printed on 06 Jun 2017**

---

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)