

## Datasheet: AHP1124

<b>Description:</b>	GOAT ANTI SNAP-25 (C-TERMINAL)
<b>Specificity:</b>	SNAP-25 (C-TERMINAL)
<b>Other names:</b>	SYNAPTOSOMAL-ASSOCIATED PROTEIN
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			2.0 - 4.0ug/ml
ELISA	▪			1/32000
Immunoprecipitation			▪	
Western Blotting	▪			0.01 - 0.03ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

**(1)Heat-mediated antigen retrieval is required prior to staining paraffin sections with this antibody. Bio-Rad recommend citrate buffer, pH6.0, for this purpose.**

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Mouse Based on sequence similarity, is expected to react with:Rat, Dog <b>N.B.</b> Antibody reactivity and working conditions may vary between species.
<b>Product Form</b>	Purified IgG - liquid
<b>Antiserum Preparation</b>	Antisera to human SNAP-25 were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
<b>Buffer Solution</b>	TRIS buffered saline
<b>Preservative</b>	0.02% Sodium Azide
<b>Stabilisers</b>	0.5% Bovine Serum Albumin
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml

<b>Immunogen</b>	Peptide sequence C-DEANQRATKMLGSG corresponding to the C-terminus of SNAP-25 (NP_003072.2, NP_570824.1).
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P60880</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">6616</a>    SNAP25    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	SNAP
<b>Specificity</b>	<p><b>Goat anti Human SNAP25 antibody</b> recognizes human synaptosomal-associated protein, 25kD (SNAP-25), also known as Super protein (Sup) or Synaptosomal-associated 25 kDa protein. SNAP-25 is a 206 amino acid, presynaptic plasma membrane protein containing 2 <a href="#">tSNARE coiled-coil</a> domains and is essential for Ca<sup>2+</sup>-regulated exocytosis of neurotransmitters. SNAP-25 binds to and negatively modulate Ca<sup>2+</sup> channels, increasing the speed of neurotransmitter release (<a href="#">Rettig et al. 1996</a>). Goat anti Human SNAP25 antibody detects an epitope within the C-terminal (CT) of human SNAP25 (aa 193-206).</p> <p>As a result of alternative splicing, SNAP-25 is expressed by neurons and neuroendocrine cells as two predominant isoforms; SNAP-25a and SNAP-25b. SNAP-25a is expressed early in development, whereas SNAP-25b expression is linked to synaptogenesis (<a href="#">Bark and Wilson 1994</a>). Goat anti human SNAP25 antibody is expected to recognise both isoforms.</p> <p>SNAP-25 may be a possible therapeutic treatment target for sufferers of hyperactivity disorders (<a href="#">Barr et al. 2000</a>), whilst SNAP-25 based peptides may be a possible treatment for involuntary muscle spasms (<a href="#">Hodel 1998</a>).</p>
<b>Western Blotting</b>	AHP1124 detects a band of approximately 22kDa in human and mouse brain lysates. A minimum incubation time of 1 hour is recommended with this antibody.
<b>References</b>	1. Hodel, A. (1998) SNAP-25. <a href="#">Int J Biochem Cell Biol. 30 (10): 1069-73.</a>
<b>Storage</b>	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Shelf Life</b>	18 months from date of despatch.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10058 available at: 10058: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

## Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025C\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M265466:150109'

**Printed on 20 Jun 2018**

---

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)