

## Datasheet: AHP1752

<b>Description:</b>	GOAT ANTI HUMAN C3
<b>Specificity:</b>	C3
<b>Other names:</b>	COMPLEMENT COMPONENT 3
<b>Format:</b>	Serum
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Turbidimetry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Serum - liquid
<b>Antiserum Preparation</b>	Antisera to human C3 were raised by repeated immunisation of goats with highly purified antigen.
<b>Buffer Solution</b>	TRIS buffered saline
<b>Preservative Stabilisers</b>	0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Immunogen</b>	Purified human C3.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P01024</a>   <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">718</a>   C3   <a href="#">Related reagents</a></p>
<b>Synonyms</b>	CPAMD1

**Specificity**      **Goat anti Human C3 antibody** recognizes human C3, a secreted protein that plays a central role in the activation of the complement system.

**Goat anti Human C3 antibody is available in bulk quantities and may be made available in other formats. Please enquire for further details.**

---

<b>Storage</b>	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Shelf Life</b>	18 months from date of despatch.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10332 available at: 10332: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10332.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10332.pdf</a>
<b>Regulatory</b>	For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M317335:180625'

**Printed on 25 Jun 2018**

---

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)